



FLEX for 6 FR Catheter



FLEX 2 mm



Notes

This drawing represents the external dimensions and the scintillator's effective point of measurement for the Medscint FLEX detector series.

For a multipoint configuration option (2CH or 4CH), all detector heads adhere to one design specifications outlined in this drawing and feature a 2-meter individual optical fiber before converging into a common cable.

For phantom insertion of the 2 mm model, Medscint recommends drilling a flat-ended cylindrical cavity with a 2.2 mm diameter, ensuring optimal detector positioning and movement.



0mm 10mm

 medscint innovative dosimetry solutions phone: 1-888-901-6521 email: info@medscint.com		
TOLERANCES UNLESS OTHERWISE SPECIFIED / SAUF INDICATION CONTRAIRE ANGLE 0 ± 1° ANGLE 0.0 ± 0.5° 0 ± 0.5mm 0.0 ± 0.1mm 0.00 ± 0.05mm		
AREA / SURFACE (mm²) N/A		
WEIGHT / POIDS (g) N/A		
TITLE / TITRE FLEX Detector Series Ext Specs		
DESCRIPTION		
SHEET / FEUILLE 1/1	PART NO. / NUMERO DE PIECE N/A	VERSION 1