


PRO / MICRO / FLASH Detector Series



Notes

This drawing represents the external dimensions and the scintillator's effective point of measurement for the Medscint detector series: PRO, MICRO, and FLASH.

For a multipoint configuration option (2CH or 4CH), all detector heads adhere to the same design specifications outlined in this drawing and feature a 2-meter individual optical fiber before converging into a common cable.

For phantom insertion, Medscint recommends drilling a flat-ended cylindrical cavity with a 5.2 mm diameter to ensure optimal detector positioning and movement.



0mm 10mm

 medscint innovative dosimetry solutions phone: 1-888-901-6521 email: info@medscint.com		
CE DOCUMENT EST LA PROPRIÉTÉ INTELLECTUELLE DE MEDSCINT INC. TOUS DROITS RÉSERVÉS. TOUTE RÉPRODUCTION, REPARTITION OU COMMUNICATION DE CE DOCUMENT, SANS LE CONSENTEMENT ÉCRIT DE MEDSCINT INC., EST STRICTEMENT INTERDITE ET POURRA ÊTRE SANCTIONNÉE.		
THIS DOCUMENT IS THE INTELLECTUAL PROPERTY OF MEDSCINT INC. ALL RIGHTS RESERVED. ANY REPRODUCTION, REPARTITION OR COMMUNICATION OF THIS DOCUMENT WITHOUT THE WRITTEN CONSENT OF MEDSCINT INC. IS STRICTLY PROHIBITED AND MAY BE SANCTIONED.		
TOLERANCES UNLESS OTHERWISE SPECIFIED / SAUF INDICATION CONTRAIRE	DIMENSIONS ARE IN MILLIMETERS	DATE: YYYY-MM-DD
ANGLE 0 ± 1° ANGLE 0.0 ± 0.5° 0 ± 0.5mm 0.0 ± 0.1mm 0.00 ± 0.05mm	DRAWN BY / DESSINE PAR B.COTE	2025-04-01
AREA / SURFACE (mm²) N/A	MATERIAL / MATÉRIEL N/A	
WEIGHT / POIDS (g) N/A	SURFACE FINISH / FINI DE SURFACE RAW	
TITLE / TITRE PRO MICRO FLASH Detector Series Ext Specs		
DESCRIPTION		
SHEET / FEUILLE 1/1	PART NO. / NUMERO DE PIECE N/A	VERSION 1

PAPER SIZE: A
8-1/2x11"